



Tuberculin Order Form

apha.tuberculinorders@apha.gov.uk

Date:

Practice Name:

Delivery Address:

Postcode:

Contact Name:

Telephone number:

Important Information

Do not mix different sizes of vials or different kits with different batch numbers at the test on a particular holding. Do not freeze tuberculin.

Item Description	Product Code	Quantity Required
Tuberculin must be used efficiently - for tests of 80 or less animals, use 2.2ml vials only		
Complete 2.2ml kit box = 400 doses in total (contains 20 vials of avian and 20 vials of bovine, 20 doses per vial)	RBI5070	
Tuberculin must be used efficiently - for tests of 81 or more animals, use 5.2ml vials only		
Complete 5.2ml kit box = 500 doses in total (contains 10 vials of avian and 10 vials of bovine, 50 doses per vial)	RBI6275	

**Approval must be given by a registered Veterinary Surgeon.
The form must either be signed below or emailed from a veterinarian.**

Approver's Name:

Approver's Signature:

Approver's RCVS Number:

Please tick this box to confirm that the approver's signature is that of a registered Veterinary Surgeon

SUBMIT

Please allow 10 working days for receipt of your order